

# THE RELIEF YOU DESERVE.

Introducing Prolieve Thermodilatation System





# The First Ever FDA Approved Microwave Thermodilatation Device for BPH is now in Singapore

Treatment of BPH just got easier! Prolieve is a first in its class microwave technology with an added balloon dilatation feature that provides immediate and lasting relief for BPH patients after just a 45-minute treatment, with a high post-treatment catheter-free rate of 85%. <sup>2</sup>

# BECAUSE WE KNOW YOU CARE

## Marketed by:



For Healthcare Professional Only. Full prescribing information will be provided upon request.

Manufactured by:



### References

Proscar) is typically indicated.

1. David Bock, et.al., 'Prolieve Transurethral Microwave Thermodilatation versus Finasteride Results of a Multicenter, Randomized Trial in Symptomatic Patients with Benign Prostatic Hyperplasia', AUA Annual Meeting, San Francisco, California, USA, November 2004.

Indication: Transurethral microwave therapy device for the treatment of symptomatic Benign Prostatic Hyperplasia (BPH) in men with a prostate size of 20-80 grams and prostatic urethral length between 1.2cm and 5.5cm and in whom drug therapy (Finasteride or

2. Bidargaddi, V & Mon, J, 'Post-Marketing Study Using Prolieve for the Treatment of Benign Prostatic Hyperplasia (BPH), 120-Month Post-Approval Study Annual Status Report, February 19, 2004, Columbia, United States.